

Patient Response to an Aerosolized Bronchodilator Administered Via the Aerogen® Solo Versus a Jet Nebulizer During Non-invasive Ventilation in Patients With Acute Exacerbation of COPD

Avdeev SN, Nuralieva GS, Soe AK, et al. Comparison of Vibrating Mesh and Jet Nebulizers During Noninvasive Ventilation in Acute Exacerbation of Chronic Obstructive Pulmonary Disease. *J Aerosol Med Pulm Drug Deliv.* 2021;34(6):358-365.

Background



There are limited data on the efficacy of aerosolized bronchodilators administered during non-invasive ventilation in patients with acute exacerbation of COPD, including the performance of different aerosol delivery devices in this setting

Objective




The aim of this study was to compare the efficacy of bronchodilator therapy delivered using the Aerogen Solo versus a jet nebulizer in patients with acute exacerbation of COPD undergoing non-invasive ventilation

Materials and Methods

Prospective randomized controlled cross-over design


Adults aged >40 years with exacerbation of COPD admitted to the respiratory unit of an acute care hospital



Non-invasive ventilation indicated because of acute hypercapnic respiratory failure*

Treatment with salbutamol 2.5 mg/2.5 mL

Randomization



*pH <7.35, partial pressure of carbon dioxide \geq 45 mmHg; [†]Operated until sputter (median 11 min); [‡]Operated until no more aerosol was produced (median 7 min). COPD, chronic obstructive pulmonary disease.


Patient Response to an Aerosolized Bronchodilator Administered Via the Aerogen® Solo Versus a Jet Nebulizer During Non-invasive Ventilation in Patients With Acute Exacerbation of COPD

Avdeev SN, Nuralieva GS, Soe AK, et al. Comparison of Vibrating Mesh and Jet Nebulizers During Noninvasive Ventilation in Acute Exacerbation of Chronic Obstructive Pulmonary Disease. J Aerosol Med Pulm Drug Deliv. 2021;34(6):358-365.




Flow and volume responses and improvements in breathlessness were superior with bronchodilator delivery via the Aerogen Solo versus a jet nebulizer in patients with acute exacerbation of COPD undergoing non-invasive ventilation

Change from baseline in FEV₁ after 120 minutes



Clinically significant* improvements in lung function with the Aerogen Solo versus a jet nebulizer



PM1107 *Refers to the between-group difference in change from baseline. BPM, breaths per minute; COPD, chronic obstructive pulmonary disease; FEV₁, forced expiratory volume in 1 second; FVC, forced vital capacity; PaCO₂, partial pressure of carbon dioxide.

Want to know more? Scan or click the QR code

